MiniLink[®] REAL-Time

TRANSMITTER


TRANSMISOR

RELEASED

A photograph of a human ear with a hearing aid device attached. The ear is shown in profile, facing right. The hearing aid is a small, light-colored device. The background is a light blue color with several geometric shapes, including a large triangle and a smaller triangle, that create a sense of depth and focus on the ear.

RELEASED

The Medtronic MiniLink® transmitter is a component of select continuous glucose sensing systems and the MiniMed® 530G insulin pump. The transmitter powers the glucose sensor, collects glucose data, and wirelessly sends the data to a MiniMed 530G insulin pump.

A complete MiniLink kit includes:

- MiniLink transmitter (MMT-7703)
- Tester (MMT-7706)
- Sensor insertion device
- Blue charger (MMT-7705)
- AAA or LR-03 alkaline battery(ies)

Indications for use

The transmitter is indicated for use as a component of select Medtronic continuous glucose sensing systems and the MiniMed 530G insulin pump.

Contraindications

None known.

Warnings

Product contains small parts and may pose a choking hazard for young children.

The sensor should be removed if redness, bleeding, pain, tenderness, irritation, or inflammation develops at the insertion site, or if you experience unexplained fever.

The optional occlusive dressing should be removed if irritation or reaction develops.

Bleeding may occur after inserting the sensor. Make sure that the site is not bleeding before connecting the transmitter to the sensor.

- If bleeding occurs, apply steady pressure with a sterile gauze or clean cloth at the insertion site until bleeding stops. After bleeding stops, connect the transmitter to the sensor.
- If bleeding persists after three minutes, remove the sensor and discard. Insert a new sensor in a different location.

Contact the 24 Hour HelpLine if you experience any adverse reactions associated with the transmitter or sensor.

Magnetic fields

Do not expose your transmitter to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. If your transmitter is inadvertently exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine for further assistance.

X-rays, MRIs, and CT scans

If you are going to have an X-ray, CT scan, MRI, or other type of exposure to radiation, disconnect your transmitter and remove your sensor before entering a room that contains any of this equipment.

Important information about airport security systems, and using your insulin pump on an airplane, can be found on the Emergency Card. Be sure to carry the Emergency Card provided when you are traveling.

Precautions

Establish a rotation schedule for choosing new sensor sites. Avoid sites that are constrained by clothing, have scar tissue, or are subject to rigorous movement during exercise.

Notice

This device complies with the United States Federal Communications Commission (FCC) and international standards for electromagnetic compatibility.

This device complies with Part 15 Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesirable operation.

The transmitter does not interfere with any radio frequency signals transmitted from outside sources. These FCC standards are designed to provide reasonable protection against excessive radio frequency interference and prevent undesirable operation of the device from unwanted electromagnetic interference.

Important: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

RF interference from other devices

Common consumer electronic devices that transmit in the same frequency band used by the transmitter may prevent the receiving device (MiniMed 530G insulin pump) from receiving the glucose information sent by the transmitter. Most cellular (mobile) phones and 900 MHz cordless phones, when transmitting or receiving, may cause significant interruption of transmitter-receiver communication. It is likely that other devices operating in similar frequency ranges will have a similar effect. This interference, however, will not cause any incorrect data to be sent and will not cause any harm to your transmitter.

MiniMed 530G insulin pumps include a programmable “Weak Signal” alert that notifies you when one or more transmitter transmissions were not received as expected by the receiving device. (The receiving device will also issue a “Lost Sensor” alert if communication is interrupted for approximately 40 minutes.)

Communication problems can typically be resolved by ensuring that the distance between transmitter and receiving device is less than six feet (1.8 meters), and by turning off or moving away from other RF transmitting devices. You can also reorient or relocate the transmitter and/or the receiving device to try to correct the interference. Testing conducted with several different cellular phones suggests that interference will not be a problem if the phone is at least 12 inches (31 cm) from the transmitter or receiving device while it is being used (greater separation distance may be required for certain devices).

Assistance

Medtronic MiniMed provides a 24 Hour HelpLine for assistance. The HelpLine is staffed with representatives who are trained in the set-up and operation of your CGM system. When calling the HelpLine, please have your pump serial number available. Your pump serial number and the 24 Hour HelpLine phone number are listed on the back of your device.

Department	Telephone number
24 Hour HelpLine (calls within the United States)	800 646 4633
24 Hour HelpLine (calls outside the United States)	+1 818 576 5555
Web site	www.medtronicdiabetes.com


Charger

The transmitter contains a non-replaceable, rechargeable battery that you can recharge as needed with the charger. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If you see a red light, see the Troubleshooting section. The charger needs one AAA alkaline battery, size E92, type LR03 (Energizer brand recommended) to operate.

Note: *A new AAA or LR-03 battery contains enough power to recharge the transmitter more than 40 times. If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.*

Installing a battery in the blue charger

- 1 Remove the charger battery cap by turning it counter-clockwise 1/4 turn using a coin in the groove of the cap.
- 2 Insert a new AAA or LR-03 battery with the flat (-) end first. Make sure that you align the small bumps on the battery cap with the small notches in the charger's battery opening. Push in the cap all the way using a coin. Turn the cap clockwise 1/4 turn to close.


Charging the transmitter


Before using the transmitter for the first time, you must fully charge the transmitter battery. This may take up to eight hours. It is recommended to recharge the transmitter after each sensor use. If you choose to recharge the MiniLink after a 6-day sensor use, the charging time will be less than 20 minutes. A fully charged transmitter battery works more than 14 days without recharging. After 14 days of use, the transmitter recharges in less than two hours.

To charge the transmitter:

- 1 If a green light on the transmitter is lit or flashing, do not connect it to the charger. The transmitter will not charge with its green light on. Wait until the green light turns off (approximately 30 seconds), and then connect the transmitter to the charger.
- 2 Connect the transmitter to the charger by lining it up, flat side down, with the charger. Push the two components together fully.

- 3 Within 10 seconds after the transmitter is connected, a green light on the charger will flash for 1-2 seconds as the charger powers on. For the rest of the charging time, the charger's green light will continue to flash in a pattern of four flashes with a pause between the four flashes.
- 4 When charging is complete, the green light on the charger will stay on, without flashing, for 15–20 seconds and then turn off.
- 5 After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter will flash for about 5 seconds and then turn off.

Setting up the transmitter

For instructions on setting up your transmitter, see the user guide for your pump.


Connecting the transmitter to the sensor

Before connecting the transmitter and sensor, make sure of the following:

- The transmitter is fully charged.
- Consult the user guide for your pump for details on setting up and starting the sensor.
- Insert the sensor and wait five minutes.
- The sensor insertion site is not bleeding.

To connect the transmitter to the sensor:


- 1 After the sensor is inserted, apply the Enlite overtape to the sensor. Consult the user guide for your Enlite senter for details on applying Enlite overtape.
- 2 Touch the rounded end of the inserted sensor to prevent it from moving during connection.
- 3 Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.
- 4 Slide the transmitter onto the sensor until the sensor's flexible arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated, the green light on the transmitter will flash within 10 seconds.

- 5 If the transmitter light does not flash, disconnect it from the sensor, wait for several seconds and then reconnect. If the transmitter light still does not flash, charge the transmitter.
- 6 After the transmitter light flashes green, use your pump to communicate with the sensor. For more instructions, see the user guide for your pump.
- 7 After the transmitter successfully sends sensor data to the pump, attach the adhesive tab from the Enlite sensor (MMT-7008) to the transmitter.
- 8 [Optional]: Apply occlusive dressing over the transmitter and the sensor.

Disconnecting the transmitter from the sensor

- 1 Carefully remove any occlusive dressing from the transmitter and sensor.
- 2 Remove the sensor's adhesive tab from the transmitter.
- 3 Hold the transmitter as shown, and pinch the flexible side arms of the sensor between your thumb and forefinger.
- 4 Gently pull the transmitter away from the sensor.

Bathing and swimming


After the transmitter and sensor are connected, they form a water-tight seal to a depth of eight feet (2.4 meters) for up to 30 minutes. You can shower and swim without removing them. No occlusive dressing is needed.

Tester

The tester is used to test the transmitter to make sure it is working.

Connecting the tester


- 1 Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter.
- 2 Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter.
- 3 Within 20 seconds the green light on the transmitter will flash for about 10 seconds when properly connected.
- 4 Check the sensor icon on the pump to ensure that the transmitter is sending a signal (see the user guide for your pump).

Disconnecting the tester


- 1 Hold the transmitter body as shown and pinch the side arms of the tester.
- 2 With the tester arms pinched, gently pull the transmitter away from the tester.

Note: To save transmitter battery life, do NOT leave the tester connected after testing.


Troubleshooting

Question: Why do I see quick flashing red lights on the charger?


Answer: The transmitter battery is very low. Leave the transmitter on the charger for 8 hours to completely recharge. If the red light is still flashing after 8 hours, leave your transmitter on the charger for 24 hours. If the red light is still flashing after its been charging continuously for 24 hours, call the HelpLine. It may be time to replace your transmitter.

Question: Why did the flashing green charger light turn off and a longer flashing red charger light turn on during charging?

Answer: Your charger battery is low. Make sure the transmitter is not connected to the charger and replace the charger battery with a new AAA or LR-03 battery.

Question: Why do I see a mix of quick and long flashing red lights on the charger?

Answer: Your charger AND transmitter batteries are very low. Replace the charger's AAA or LR-03 battery. If you now get the pattern for very low transmitter battery, leave the transmitter on the charger for 8 hours to recharge. If the red light is still flashing after 8 hours, leave your transmitter on the charger for 24 hours. If the red light is still flashing after its been charging continuously for 24 hours, call the HelpLine. It may be time to replace your transmitter.

Question: I had my transmitter on the charger for a day. Will this damage my transmitter?

Answer: It will not damage the transmitter. You cannot overcharge it.

Question: What should I do if the transmitter green light does not flash when connected to the sensor?

Answer: Is the sensor inserted in the body? **If it is not inserted**, the transmitter will not flash green or send signals to an insulin pump.

If the sensor is inserted in the body, you need to disconnect the transmitter from the sensor, wait for several seconds and then reconnect. If the green light still does not flash, charge the transmitter.

Question: Why didn't I see the transmitter green light flash after connecting it to the tester?

Answer: Check the connection. If you still do not see a green light flash, fully recharge the transmitter battery. Test the transmitter with the tester. If you still do not see a green light flash, call the HelpLine. It may be time to replace your transmitter.

Cleaning the transmitter

Caution: **Do not discard the transmitter in a medical waste container or otherwise subject it to incineration. The transmitter contains a battery that may explode upon incineration.**

Caution: **The charger and tester are NOT watertight. Do NOT immerse in water or any other cleaning agent. If liquid gets inside the black connector of the tester, shake the liquid out and allow to air dry.**

For single-patient use

- 1 Wash your hands thoroughly.
- 2 Attach the tester to the transmitter.
- 3 Dampen a clean cloth with mild liquid soap and warm water. Wipe the outside of the transmitter.
- 4 Rinse the transmitter under warm tap water.

- 5 Using an antibacterial hand-sanitizer (available at a local drugstore) on a clean, dry cloth, wipe the transmitter's surface. Do NOT get any hand-sanitizer inside the tester opening or transmitter connector. Repeated exposure to hand-sanitizer could damage the connectors and affect the transmitter's performance as a result. If you get hand-sanitizer inside the tester or connector, allow them to air dry.
- 6 Disconnect the tester from the transmitter.
- 7 Place the transmitter on a clean, dry cloth and air dry for 2–3 minutes.

For multiple-patient use

Caution: The charger and tester are NOT watertight. Do NOT immerse in water or any other cleaning agent. If liquid gets inside the black connector of the tester, shake the liquid out and allow to air dry.

Warning: Always clean the transmitter after removing it from the patient and before attaching it to the charger. If there is any blood inside the connector, the transmitter must be discarded.

Dispose the transmitter according to the local regulations for battery disposal (non- incineration).

Follow this procedure to clean the transmitter after each patient use:

- 1 Wash your hands thoroughly.
- 2 Attach the tester to the transmitter.
- 3 Dampen a clean cloth with mild liquid soap and warm water. Wipe the outside of the transmitter.
- 4 Rinse the transmitter under warm tap water.
- 5 Apply 3-4 drops of a quaternary ammonium compound disinfectant (for example, Cavicide®) on a clean, dry cloth and wipe the transmitter.
- 6 Hold the transmitter and wipe with 70% Isopropyl alcohol.
- 7 Disconnect the tester from the transmitter.
- 8 Place the transmitter on a clean, dry, non-shedding cloth and air dry.

- 9 Once it is completely dry, place it in a sealed bag labeled with the cleaning date (see the following label example).

Device:
Date:
Method of Decontamination:
Wash:
Disinfectant:

Cleaning the charger or tester

Caution: The blue charger and tester are not water-tight. Do NOT immerse in water or any other cleaning agent. Do not allow liquid to come in contact with the blue charger connector. Repeated exposure to liquid could damage the connector and affect the performance of the charger. If liquid comes in contact with the connector, allow the charger to air dry before proceeding with the cleaning.

For single-patient use

- 1 Wash your hands thoroughly.
- 2 Use a damp cloth with mild cleaning solution, such as Dial® Dishwashing detergent, to clean any dirt or foreign material from the outside of the charger or tester. Never use organic solvents such as paint thinner or acetone to clean the charger or tester.
- 3 Using an antibacterial hand sanitizer (available at a local drugstore) on a clean, dry cloth, wipe the outside surface areas of the charger or tester. Do NOT get any hand-sanitizer inside the tester opening or on the black connector. If you get hand sanitizer inside these locations, allow it to air-dry.
- 4 Place the charger or tester on a clean, dry cloth and air dry for 2-3 minutes.

For multiple-patient use

Caution: Dispose the charger according to the local regulations for battery disposal (non-incineration).

Warning: Always clean the transmitter after removing it from the patient and before attaching it to the charger. If blood comes in contact with any surface areas of the charger or tester, the contaminated device must be discarded. The charger contains a battery which may explode upon incineration.

Follow this procedure to clean the charger or tester after each patient use:

- 1 Wash your hands thoroughly.
- 2 Use a damp cloth with mild cleaning solution, such as Dial Dishwashing detergent, to clean any dirt or foreign material from the outside of the charger or tester. Never use organic solvents such as paint thinner or acetone to clean the charger or tester.
- 3 Apply three to four drops of a quaternary ammonium compound disinfectant (for example, Cavicide) on a clean, dry cloth and wipe all the surface areas of the charger or tester.
- 4 Slightly dampen a clean cloth with isopropyl alcohol and wipe all the surface areas of the charger or tester.
- 5 Place the charger or tester on a clean dry, non-shedding cloth and allow to air dry.
- 6 Once it is completely dry, place it in a sealed bag labeled with the cleaning date (see the following label example).

Device: Date: Method of Decontamination: Wash: Disinfectant:
--

Storing the devices

Store the transmitter, charger, and tester in a clean, dry location at room temperature. Although not required, you may store the transmitter on the charger. If the transmitter is not in use, you must charge the transmitter at least once every 60 days.

Specifications

Biocompatibility	Transmitter: Complies with ISO 10993-1 for body contact
-------------------------	---

Applied parts	Transmitter Sensor
Operating conditions	Transmitter Temperature: +32° to +122°F (0° to 50°C) Caution: When operating the transmitter on a tester in air temperatures greater than 106°F (41°C), the temperature of the transmitter may exceed 109°F (43°C). Transmitter Relative Humidity: 10% to 95% with no condensation Charger Temperature: +50° to +104°F (10° to 40°C) Charger Relative Humidity: 30% to 75% with no condensation
Storage conditions	Transmitter Temperature: -4° to +131°F (-20° to +55°C) Transmitter Relative Humidity: 10% to 100% with condensation Charger Temperature: +14° to +122°F (-10° to +50°C) Charger Relative Humidity: 10% to 95% with no condensation
Battery life	Transmitter: 14 days of continuous glucose monitoring immediately following a full charge Charger: Completes 40 typical charging operations using a new AAA or LR-03 battery
Transmitter frequency	MMT-7703NA 916.5 megahertz MMT-7703WW 868.35 megahertz

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The MiniLink transmitter (MMT-7703 transmitter) is intended for use in the electromagnetic environment specified below. The customer or the user of the MiniLink transmitter should make sure that it is used in such an environment. The transmitter contains no repairable parts.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The MiniLink transmitter uses RF energy only for system communication functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MiniLink transmitter is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	<p>Note: <i>The preceding statement is required by IEC 60601-1-2 for Group 1, Class B devices. However, since the MiniLink transmitter is battery powered, its emissions will not be affected by the establishment power supply and there is no evidence of any issues associated with the use of the system in domestic establishments.</i></p>
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV indirect ±8 kV air	±6kV indirect (30-60% relative humidity) ±30 kV air (<5% relative humidity)	The MiniLink transmitter should not be affected by electrostatic discharge that might occur under normal conditions of use.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for input/output lines	Not applicable Not applicable	
Surge IEC 61000-4-5	±1 kV lines(s) to line(s) ±2 kV line(s) to earth	Not applicable Not applicable	
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	Not applicable Not applicable Not applicable Not applicable	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the MiniLink transmitter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> <p>$d=1.2\sqrt{P}$ 80MHz to 800 MHz $d=2.3\sqrt{P}$ 800MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>NOTE: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) tele-phones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MiniLink transmitter is used exceeds the applicable RF compliance level above, the MiniLink transmitter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the MiniLink transmitter.</p>			

RELEASED

Recommended separation distances between portable and mobile RF communications equipment and the MiniLink transmitter (MMT 7703 transmitter)		
This section provides information on the recommended separation distance between portable and mobile RF communications equipment and the MiniLink transmitter. The MiniLink transmitter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MiniLink transmitter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MiniLink transmitter as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter (W)	Separation distance according to the frequency of transmitter (m)	
	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.5GHz $d=2.3\sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.74
1	1.2	2.3
10	3.8	7.4
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Warranty

Medtronic MiniMed warrants the Medtronic MiniLink transmitter to the purchaser of the product against defects in material and workmanship for a period of twelve (12) months and the charger for up to six (6) months from the date of purchase.

During the warranty period, Medtronic MiniMed will repair or replace, at its discretion, any defective MiniLink transmitter or charger, subject to the conditions and exclusions stated herein. This warranty applies only to new devices. In the event a MiniLink transmitter or charger is repaired or replaced, the warranty period will not be extended past its original expiration date.

This warranty is valid only if the Medtronic MiniLink transmitter or charger is used in accordance with the manufacturer's instructions. Without limitation, this warranty will not apply:

- If damage results from changes or modifications made to the MiniLink transmitter or charger by the user, or third parties, after the date of sale;


- If service or repairs are performed by any person or entity other than the manufacturer;
- If damage results from a *Force Majeur* or other event beyond the control of the manufacturer;
- If damage results from negligence or improper use, including but not limited to: improper storage, submersion in fluid, physical abuse (such as dropping); or
- If blood or water has entered the inside of the MiniLink transmitter connector.

This warranty shall be personal to the original user. Any sale, rental or other transfer or use of the product covered by this warranty to or by a user other than the original user shall cause this warranty to immediately terminate. This warranty does not apply to Glucose Sensors and other accessories.

The remedies provided for in this warranty are the exclusive remedies available for any defects in material or workmanship in the product. Any statutory rights granted to consumers under any applicable legislation are reserved. Neither Medtronic MiniMed nor its suppliers or distributors shall be liable for any incidental, consequential, punitive or special damages of any nature or kind caused by or arising out of a defect in the product. All other warranties, except any applicable mandatory statutory warranties, expressed or implied, are excluded and specifically disclaimed, including, but not limited to, any warranty of merchantability or fitness for a particular purpose.

RELEASED

Icon Table

SN	Serial number
REF	Catalogue number
(1X)	One per container/package
	Date of manufacture (year - month)
	Manufacturer
	Follow instructions for use
	Storage temperature range
	Radio communication device
CONF	Configuration
	Type BF equipment (protection from electrical shock)
IPX8	Transmitter: Protected against the effects of continuous immersion in water (8 feet or 2.4 meters immersion for 30 minutes).
	Attention: Read all warnings and precautions in instructions for use.
	Storage humidity range

©2012 Medtronic MiniMed, Inc. All rights reserved.

MiniMed® and MiniLink® are registered trademarks of Medtronic MiniMed, Inc. Enlite™ is a trademark of Medtronic MiniMed, Inc. Cavicide® is a registered trademark of Metrex. Dial® is a registered trademark of Henkel Corporation.

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©2012 Medtronic MiniMed, Inc. Reservados todos los derechos.