**Medtronic**

ENLITE™

Glucose Sensor  
User Guide

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
RELEASED

## Introduction

The Enlite™ glucose sensor (sensor) is part of your Continuous Glucose Monitoring System. The sensor continuously converts tiny amounts of glucose from the fatty layer under your skin into an electronic signal. Your system then uses these signals to provide sensor glucose values.

- A. Pedestal
- B. Needle housing
- C. Sensor
- D. Clear liner

## Potential risks

### Risks related to sensor use

General risks with sensor use include:

- Skin irritation or reaction to adhesives
- Bruising
- Discomfort
- Redness
- Bleeding
- Pain
- Rash
- Infection
- Irritation from tapes used with glucose-sensing products.
- Raised bump
- Appearance of a small “freckle-like” dot where needle was inserted
- Allergic reaction
- Fainting secondary to needle insertion
- Soreness or tenderness
- Swelling at insertion site
- Sensor fracture, breakage or damage
- Minimal blood splatter associated with sensor needle removal
- Residual redness associated with adhesive and or tapes
- Scarring

### Specific risks related to sensor use

During the conduct of the Performance Evaluation of the Enlite™ Glucose Sensor to Support a Full 144 Hours (6 Days) of Use<sup>1</sup>, the following specific sensor risks were identified:

<sup>1</sup> Medtronic Inc., A Performance Evaluation of the Enlite™ Glucose Sensor to Support a Full 144 Hours (6 Days) of Use, CER247/Z25/C, May 15, 2012.

- Pain at sensor insertion site during sensor wear (1 incident reported)
- Dermatological risks (number of events/number of skin assessments performed):
  - Redness due to device insertion (abdomen, 35/213)
  - Redness in area of adhesive (abdomen, 63/213)
  - Bruising (abdomen, 2/213)
  - Bump (induration) (abdomen, 1/213)
  - Bleeding:
    - Bleeding immediately after insertion (abdomen, 2/213)
    - Bleeding at removal (abdomen, 5/213)
  - Fluid discharge (abdomen, 1/213)
  - Other, including indentation from device, for example, the transmitter (abdomen, 9/213)

When using medications containing acetaminophen some people experience a temporary increase in sensor glucose readings for less than one hour. For some people the sensor readings are not noticeably affected by acetaminophen. Always use BG meter readings to verify your glucose level before making therapy decisions.

### Indications for use

The sensor is intended for use with Medtronic Diabetes (Medtronic) glucose sensing systems to continuously monitor glucose levels in persons with diabetes.

### Contraindications

None known.

### General warnings

Do not make therapy treatment decisions based on sensor glucose values because sensor glucose and blood glucose values may differ. If your sensor glucose reading is high or low, or if you feel symptoms of high or low glucose, confirm your blood glucose with your BG meter prior to making therapy decisions to avoid severe high or low glucose conditions.

Keep an arm's length between you and the patient when removing the needle. A retractable needle is attached to the sensor and minimal blood splatter may occur.

Keep the needle housing within sight at all times to avoid accidental needle stick injury.

Always inspect the packaging for damage prior to use. Sensors

are sterile and nonpyrogenic, unless the package has been opened or damaged. Do not use any of the sensors if the sterile package has been opened or damaged. Use of an unsterile sensor can cause site infection.

Watch for bleeding at the insertion site (under, around, or on top of the sensor). Apply steady pressure using sterile gauze or a clean cloth placed on top of the sensor for up to three minutes, if bleeding occurs. The use of unsterile gauze can cause an infection.

Do not attach the transmitter to the sensor, if bleeding does not stop. Blood can get into the transmitter's connector, which can damage the device. Remove the sensor and continue to apply steady pressure until the bleeding stops.

1. Remove the sensor and discard.
2. Check the site for redness, bleeding, irritation, pain, tenderness, or inflammation and treat accordingly.
3. Insert a new sensor in a different location.

Keep the sensor away from children. This product contains small parts and may pose a choking hazard.

### General precautions

Wash your hands with soap and water before inserting the sensor to help prevent site infection.

Wear gloves when inserting the sensor into someone other than yourself to avoid contact with patient blood.

Do not insert the sensor through tape to ensure proper sensor function and insertion.

Only use alcohol to prepare the insertion site to ensure that residue is not left on the skin.

Rotate the sensor insertion site so that sites do not become overused.

Discard used sensors and needle housings in a sharps container after each use to avoid accidental needle stick injury.

Do not clean, resterilize, or try to extract the needle from the needle housing. An accidental needle stick may occur.


Do not reuse sensors. Reuse of the sensor may cause damage to the sensor surface and lead to inaccurate glucose values, site irritation, or infection.

## Where to insert the sensor

**WARNING:** Always refer to the Enlite Serter User Guide for instructions on inserting the sensor. Failure to follow directions may result in pain or injury.

**CAUTION:** Avoid the 2 inch (5.0 cm) area around the navel to help ensure a comfortable insertion site, and to help with adhesion.

Choose an insertion site that has an adequate amount of subcutaneous fat. Shown here are the best body areas (shaded) for sensor insertion.

Do not insert the sensor in muscle or areas constrained by clothing or accessories; areas with tough skin or scar tissue; sites subjected to rigorous movement during exercise; or in sites under a belt or on the waistline to avoid chafing or site irritation.

## Removing the sensor

When you are ready to change your sensor, disconnect the transmitter from the sensor as described in your transmitter user guide. Gently pull the sensor from your body to remove it. Place it in a sharps container.

## Reagents

The sensor contains two biological reagents: glucose oxidase and human serum albumin (HSA). Glucose oxidase is derived from *Aspergillus Niger*, purified and dried according to Type VII-S guidelines and cross-linked to HSA, purified and dried albumin fraction V, derived from pasteurized human serum with glutaraldehyde. Since less than 0.4  $\mu\text{g}$  of glucose oxidase and less than 0.7  $\mu\text{g}$  of HSA is used to manufacture each sensor, the risks of tissue reactions and viral transmission are considered minimal.

## Storage and handling

**CAUTION:** Do not freeze the sensor, or store it in direct sunlight, extreme temperatures, or humidity. These conditions may damage the sensor.

Store sensors at room temperature between +36°F to +86°F (+2°C to +30°C). If the sensor is stored in a cool environment, allow the sensor to warm to room temperature to prevent condensation. Do not store sensors at temperatures lower than +36°F (+2°C).

## MP6025676-2AF1 / A

Discard sensor after “Use by” on label, or if the package is damaged, or the seal is broken.

### Specifications

Dimensions	2.46 x 2.50 x 1.87 inches
Weight	6.25 grams

### Sensor life of use

The sensor can be used one time, and it has a maximum life of 144 hours (six days). The 144-hour life span of the sensor begins when the insulin pump or monitor receives the first METER BG NOW alert.

### Icon table

	Use by
	Do not reuse
	Follow instructions for use
REF	Catalogue number
	Batch code
	Sterilized using irradiation
	Temperature limitation
	Date of manufacture
	Open here
	Manufacturer

©2012 Medtronic MiniMed, Inc. All rights reserved. Enlite™ is a trademark of Medtronic MiniMed, Inc.

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